# ProCheck

---

**User Manual**

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# Intended Use

EN

ProCheck is designed for health data management and monitoring via mobile devices. It provides an easy-to-use interface to track Blood Glucose with you MediSure Empower Blood Glucose Monitor.

Not only logging your daily measurement data, ProCheck also provide analytic function to gather your testing records, and converts them to visualized graphs and charts. This application provides organized self-monitoring information to you and your healthcare provider.

V3.0 2022/06  
541-0099100-009

# Home

When you launch the application, you will see the latest blood testing result with date and time.

You can check testing result of blood glucose by tapping blood glucose button.

## Home

### Blood Glucose Data

- Display most recent data and time for blood glucose
- Background color will change

Red: above High Target  
Yellow: below Low Target  
(‘Setting’ – ‘Health Target’)

### Measurement Button

- Start headphone plug Dongle measurement

### Blood Glucose Button

- Tap to display most recent data for blood glucose

### Blood Glucose Data List

- Display most recent data for blood glucose

### TabBar

- Tap to switch to 'Home', 'Data', 'Analysis', 'Setting', Page

### Target

- Display target range set in 'Setting'-'Telehealth Target' page

Rotate phone in anticlockwise direction to Home

Rotate phone in clockwise direction to Summary page

Home page

Glucose				mmol/L
	GEN	AC	PC	
Target Range	3.9-5.6	3.9-5.6	3.9-7.8	
7 Days	6.4	4.9	7.1	
14 Days	5.7	5.2	7.3	
30 Days	5.7	5.8	7.8	
60 Days	5.7	6.0	8.0	

GEN: not specified  
AC: before meal  
PC: after meal

### Period

- Specify day count from now, that data during the period are included for average

### Average

- Display average of glucose values in specific period, for AC, PC, GEN, respectively
- QC data will not be included

# Data

The Data screen shows your existing test results. You can quickly review every result's information, which is arranged by date from the current date on top to previous dates while going down to the next test result.

## Period Tab

- Tap to select period want to display
- Period include 7Day, 14Day, 1 Month, 2 Month 3Month

## Data Table

- Data of blood glucose during selected period are displayed
- The latest data are displayed on top
- Display data contents include date/time, data type, value, unit, upload and note status
- data background color will change

Red: above High target

Yellow: below Low target

('Setting' - 'Health Target')

## Blood Glucose Tab

- Tap to display data of blood glucose

## Upload and Note Status


- Display Upload Status  
Uploaded   
Not Uploaded 
- Display Note status with Note   
Without Note 

Tap each record to review details or add note

## Note

- Click to add or modify note

# Analysis

ProCheck provides a lot of statistics in form of charts and graphs. You can view your collected data visually and understand it better. To check your averaging test results, you can tap each record. It takes you to the screen which shows you some personal stats and the statistics for the last 7 days/14 days/1 month/2 months/ 3 months.

## Analysis

### Period Tab

- Tap to select period want to display
- Period include 7Day, 14Day, 1 Month, 2 Month 3Month

### Blood Glucose Display

- Display blood glucose

### Average Value Display

- Average value for blood glucose during selected period
- For glucose, Before meal (AC), After meal (PC), Not specified (GEN) data is calculated separately

### Graph Button

- Click to display trend graph for blood glucose during selected period

### Blood Glucose Tab

- Tap to display data of blood glucose

### High/Normal/Low Percent

- Percentage for blood glucose during period
  - higher than target range
  - within target range
  - lower than target range

### Highest and Lowest Value

- Highest and Lowest value for blood glucose during selected period
  - Red: above High Target
  - Yellow: below Low Target

# Analysis – Graph

**Blood Glucose Information**

- Blood Glucose Unit
- For glucose, the representative symbol for AC, PC GEN

**Scrollable Graph**

- Scroll to Right: display data in newer period
- Scroll to Left: display data in older period

**Data Point and Line**

- Display Individual point and connected line for the blood glucose during

GEN are displayed separately

# Setting

## Bluetooth Smart Meter

### Bluetooth Pair Rules

- Search for Bluetooth meters automatically when entering Bluetooth Smart Meter page

### Back Button

- Go back to setting root page

### Paired Bluetooth Meter List

- List of Bluetooth meters finished pairing with App

### Scan

- Pairing BLE meter by scanning QR Code
- Tap to turn on the camera and zoom-in
- Scan meter QR Code to get Bluetooth MAC address of this meter
- Show this Bluetooth MAC address on the list of Current Meters
- Save and leave to complete pairing process
- When this meter connects to the app, replace the Bluetooth MAC address with Bluetooth Device Name

### Search

- Search for Bluetooth meters to pair
- If there is Bluetooth meter advertising in range, then pop-up Available Meters window with available meters to be paired
- Tap Add to pair this meter
- Save and leave to complete pairing process

### Paired Bluetooth Meter List

- List of Bluetooth meters finished pairing with App

## Export

### Back Button

- Go back to setting root page

### Blood Glucose and Period

- Choose period to be exported

### Email Address

- Input Email address where the exported data is sent to

### Export Data Format

- Choose export data format  
JPG – screen shot for pages  
CSV

### Export Button

- Click to export data

## TeleHealth Service

### Back Button


- Go back to setting root page

### Upload Button

- Click the button to upload data to the selected service URL site

### Last Upload Status

- Display the date and time for last upload
- No data ever upload display 'No upload records'

### Telehealth Service URL

- Tap to open TeleHealth Service list
- If the server to be chosen is not on the list, then select Others to key-in server url
- Default setting: Blank

### Automatic Upload Switch

- User could turn on/off switch, If switch is turned on, App will upload data to the selected service URL site automatically when new measurement data is imported to App
- App will not retry upload if upload attempt failed, and will run next auto-upload when new data imported
- Default setting: off

## Profile

### Back Button

- Go back to setting root page

**Profile**

**Gender**

**Birthday**

20 - 1 - 1978

**Height**

170 cm Metric

150 160 170 180 190

Home Data Analysis Setting

### Gender

- Select user gender

### Birthday

- Input user birthday

### Height

- Input user height
- Two units could be selected: Metric or Imperial

## Health Target

### Back Button

- Go back to setting root page

Glucose					
AC		3.9	-	5.6	mmol/L
PC		3.9	-	7.8	mmol/L
GEN		3.9	-	5.6	mmol/L

### Glucose Target


- Units: fixed mmol/L
- Low target and High target are specified individually for AC, PC, GEN

## About

### Back Button

- Go back to setting root page

Copyright and Software Version Information

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The screenshot shows a mobile application interface titled "Health Target". At the top left, there is a back button icon (a hand pointing left) with a blue arrow pointing to it from the text "Bouton Retour". Below the title, the section is labeled "Glucose". It contains a table with three rows of data:

Category	Icon	Lower Target (mmol/L)	Upper Target (mmol/L)	Unit
AC	Apple icon	3.9	5.6	mmol/L
PC	Person icon	3.9	7.8	mmol/L
GEN	Group of people icon	3.9	5.6	mmol/L

At the bottom of the screen is a navigation bar with four icons: Home (house), Data (bar chart), Analysis (pie chart), and Setting (gears). The Setting icon is highlighted in blue. A blue arrow points from the text "Cible glycémie" to the "mmol/L" unit in the table.

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
**[www.MediSure/support](http://www.MediSure/support)**

**MediSure Canada Inc.**

785 Westney Rd. S Unit 23,  
Ajax, ON L1S 7G1 CANADA

[www.medisure.ca](http://www.medisure.ca)

Toll Free : 1-855-634-7873

(8:30 am-5:00 pm EST, Monday-Friday)

For assistance outside of these hours, please contact your healthcare professional.

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